



Cordless MOMICARE is convenient for massaging while doing other things in any situation!

- Massaging while working
- Massaging during household tasks
- Massaging while commuting or long-distance traveling by train
- Massaging while relaxing at home

maxell  
Within, the Future

2018.3  
Maxell MXTS-MR200 Low-frequency  
Pulse Massager Catalog

Cordless Low-frequency Pulse Massager

momicare



Just place it on  
your skin!  
Easy massage

When you're uncomfortable, just slap it on.

New "MOMICARE",  
the evolution of  
convenient massage.



● Connection cable (sold separately)



White

MXTS-200CBL1P  
JAN code: 4902580756666

Major specifications

Product name	Connection cable
Model	MXTS-200CBL
Set contents	Connection cable x 1 User manual
Cord length	About 1.2 m
Weight	About 14 g



Safety precautions

To use this product securely and safely,  
be sure to read the User Manual thoroughly before use.



Danger

Never use this product in combination with any of the electric medical equipment described below.

1) Self-contained medical equipment such as a pacemaker, 2) life-support electric medical equipment such as an artificial heart-lung machine, and 3) wearable electric medical equipment such as an electrocardiograph.  
Failure to follow the above precautions may result in an accident or significant bodily harm.



Warnings

If any of the following conditions apply to you, consult a doctor before using this product.

1) Acute disease, 2) malignant tumor, 3) infectious disease, 4) pregnancy, menstruation, or immediately after having given birth, 5) cardiac disorder, 6) body temperature of 38°C or higher (fever) [during acute onset of symptoms (fatigue, chill, blood pressure change, etc.) and debilitating condition], 7) slight anemia, 8) abnormality (infection, wound, etc.) in the skin of the area to be massaged, 9) abnormal blood pressure, 10) when requiring rest or when feeling very ill, 11) acute [pain] disorder such as spine fracture, ligament rupture, or torn muscle, 12) perceptual disorder due to severe peripheral circulatory failure such as diabetes, 13) particular feelings of physical discomfort other than the above.

Failure to do so may result in an accident or poor health.

- Be sure to turn off the power when removing the gel sheet. Failure to do so may result in a strong electric shock.
- If you feel unwell while using the product, stop using it immediately and consult a doctor. Failure to do so may result in an accident or poor health.

- The product shown in this catalog complies with the Restriction of Hazardous Substances (RoHS) in EU.
- The specifications and design of the product shown in this catalog are subject to change without notice due to modification.
- The actual product may differ slightly from the photographs shown in this catalog, due to printing requirements.
- The names, logos, and service marks contained in this catalog are trademarks or registered trademarks of Maxell and other companies.
- The information contained in this catalog is accurate as of March 2018.

Contacts

<http://www.maxell.co.jp>

Maxell Contact Us

Search

■ Maxell product dealer



For global environmental protection, this catalog uses vegetable-oil-based ink.

Effectiveness/  
Efficacy

Relieves stiff shoulders, prevents paralyzed muscle from atrophying,  
and provides effective massages. This product is for home use only.

Home-use low-frequency massage device

Controlled medical device