

Cordless MOMICARE is convenient for massaging while doing other things in any situation!


- Massaging while working
- Massaging during household tasks
- Massaging while commuting or long-distance traveling by train
- Massaging while relaxing at home


○ Connection cable (sold separately)




White
MXTS-200CBL1P
JAN code: 4902580756666

Major specifications	
Product name	Connection cable
Model	MXTS-200CBL
Set contents	Connection cable x 1 User manual
Cord length	About 1.2 m
Weight	About 14 g

 Safety precautions	To use this product securely and safely, be sure to read the User Manual thoroughly before use.
---	---

 **Danger** Never use this product in combination with any of the electric medical equipment described below.

1) Self-contained medical equipment such as a pacemaker, 2) life-support electric medical equipment such as an artificial heart-lung machine, and 3) wearable electric medical equipment such as an electrocardiograph.
Failure to follow the above precautions may result in an accident or significant bodily harm.

 **Warnings** If any of the following conditions apply to you, consult a doctor before using this product.

- 1) Acute disease, 2) malignant tumor, 3) infectious disease, 4) pregnancy, menstruation, or immediately after having given birth, 5) cardiac disorder, 6) body temperature of 38°C or higher (fever) [during acute onset of symptoms (fatigue, chill, blood pressure change, etc.) and debilitating condition], 7) slight anemia, 8) abnormality (infection, wound, etc.) in the skin of the area to be massaged, 9) abnormal blood pressure, 10) when requiring rest or when feeling very ill, 11) acute [pain] disorder such as spine fracture, ligament rupture, or torn muscle, 12) perceptual disorder due to severe peripheral circulatory failure such as diabetes, 13) particular feelings of physical discomfort other than the above.
Failure to do so may result in an accident or poor health.
- Be sure to turn off the power when removing the gel sheet. Failure to do so may result in a strong electric shock.
 - If you feel unwell while using the product, stop using it immediately and consult a doctor. Failure to do so may result in an accident or poor health.

- The product shown in this catalog complies with the Restriction of Hazardous Substances (RoHS) in EU.
- The specifications and design of the product shown in this catalog are subject to change without notice due to modification.
- The actual product may differ slightly from the photographs shown in this catalog, due to printing requirements.
- The names, logos, and service marks contained in this catalog are trademarks or registered trademarks of Maxell and other companies.
- The information contained in this catalog is accurate as of March 2018.

Contacts <http://www.maxell.co.jp>


■ Maxell product dealer



maxell
Within, the Future

2018.3
Maxell MXTS-MR200 Low-frequency
Pulse Massager Catalog

Cordless Low-frequency Pulse Massager



Just place it on
your skin!
Easy massage

When you're uncomfortable, just slap it on.

New "MOMICARE",
the evolution of
convenient massage.



[Effectiveness/
Efficacy] Relieves stiff shoulders, prevents paralyzed muscle from atrophying,
and provides effective massages. This product is for home use only.

Home-use low-frequency massage device Controlled medical device

MOMICARE includes 5 modes and 10 strength settings. MOMICARE eases stiffness, pain, and fatigue.



Simple and easy to use

mode_1 Automatic

A 15-minute massage course consisting mainly of “Knead” and “Press” in addition to “Tap” and “Vibrate”.

mode_2 Quick

A 7-minute massage course consisting of “Tap”, “Vibrate”, “Knead”, and “Press”. This course is appropriate for quickly easing fatigue.

Apply MOMICARE wherever it is needed

mode_3 Neck and shoulders

A massage course consisting of “Tap” and “Vibrate” to ease your neck and shoulders. This course also eases upper arm fatigue.



mode_4 Waist

A massage course consisting of “Tap”, “Vibrate”, and “Knead” to ease lassitude, heaviness, and fatigue of waist.

mode_5 Legs

A massage course consisting of “Knead”, “Press”, and “Vibrate” to ease thigh and calf fatigue.

MOMICARE employs a frequency range of 1 to 100 Hz that provides highly effective treatment and minimal muscle fatigue.



MOMICARE employs a soft cushion gel sheet.

The cushion gel sheet has a smooth feel that makes it pleasant and enjoyable to use. It can be removed and replaced, and multiple sheets can be used – one for each family member.

* Use the MXTS-200GEL series gel sheets as replacement sheets.



MOMICARE gives priority to ease of use and safety.

MOMICARE’s thoughtful design enables you to turn it on and set the mode and strength by hand before attaching it where it is needed. The body detection sensor automatically stops MOMICARE as soon as the gel sheet leaves your skin.

Dedicated gel sheets* for each part of your body.

We recommend using these when you want to intensively massage the soles of your feet or massage a wide area such as your waist.

* Spot type/Wide type



Cordless and rechargeable

MOMICARE is a compact and cordless massage device with a built-in lithium-ion battery. A recharging power supply adapter and USB cable are included. Once the battery is charged, you can use MOMICARE 7 times for about 15 minutes*.

The connection cable (sold separately) enables you to control MOMICARE by hand.

* When the strength level is set to 10 (maximum). The number of repetitions varies according to mode.



Product range

Main unit

White

1-unit pack | MXTS-MR200W1P JAN code: 4902580756550

2-unit pack | MXTS-MR200W2P JAN code: 4902580756581

Pink

1-unit pack | MXTS-MR200P1P JAN code: 4902580756567

2-unit pack | MXTS-MR200P2P JAN code: 4902580756598

Black

1-unit pack | MXTS-MR200B1P JAN code: 4902580756574

2-unit pack | MXTS-MR200B2 JAN code: 4902580756604

Major specifications			
Brand name	MXTS Low-frequency Pulse Massager	Set contents (1-unit pack)	Main unit x 1
Model	MXTS-MR200		Gel sheet x 1
Category	Equipment 78 Home-use Electric Massage Device		Protective sheet x 1
Medical device category	Controlled medical device		Power supply adapter x 1
General name	Home-use low-frequency massage device		USB cable x 1
Medical device ID number	229AKBZX00042000	Set contents (2-unit pack)	User manual
Rated voltage of main unit	3.7 V DC (lithium-ion battery)		Main unit x 2
Dedicated power supply adapter (for charging)	Input: 100 V AC, 50 to 60 Hz, 0.4 W		Gel sheet x 2
Maximum output current	Output: 5 V DC, 0.04 A		Protective sheet x 1
Maximum output voltage	10 mA or less		Power supply adapter x 1
Power consumption (while working)	About 42 V at maximum setting	Material	USB cable x 1
Power consumption (while charging)	About 0.1 W		User manual
Basic frequency	About 0.2 W		Main unit case: PC + ABS
Rated time	About 1 to 100 Hz		Gel sheet: Acrylic gel
Output pulse width	About 15 minutes		Main unit: 44 (height) x 44 (width) x 10.5 (depth) mm
Timer	200 μ seconds	Dimensions (1 unit)	Gel sheet: 61 (height) x 125 (width) x 0.5 (thickness) mm
	Automatically turns off after about 15 minutes.		Weight (1 unit)
Charging time	Quick mode: Automatically turns off after about 7 minutes.		About 25 g (main unit + gel sheet)
3 hours (at 10 to 35°C room temperature)			

Replacement/dedicated gel sheets (sold separately)

① Standard

② Spot type

③ Wide type

① MXTS-200GELW2P
MXTS-200GELP2P
MXTS-200GELB2P

JAN code: 4902580756611
JAN code: 4902580756628
JAN code: 4902580756635

② MXTS-200GELFW2P
③ MXTS-200GELWW1P

JAN code: 4902580756659
JAN code: 4902580756642

Major specifications		Standard	Spot type	Wide type
Brand name	MXTS Gel sheet			
Model	MXTS-200GEL			
Medical device category	General medical device			
General name	Inductor for home-use low-frequency pulse massager			
Medical device registration number	40B2X10023000001			
Color	White/Pink/Black	White	White	White
Set contents	Gel sheet x 2	Gel sheet x 2	Gel sheet x 1	
	Protective sheet x 1			
	User manual			
Material	Acrylic gel			
Dimensions	61 (height) x 125 (width) x 0.5 (thickness) mm	59 (height) x 190 (width) x 0.5 (thickness) mm	154 (height) x 180 (width) x 0.5 (thickness) mm	
Weight	About 12 g (1 piece)	About 12 g (1 piece)	About 41 g (1 piece)	